



1.2 锆石 LA-ICP-MS 微区原位 U-Pb 定年和微量元素分析 (91500)

锆石 U-Pb 同位素定年和微量元素含量在武汉上谱分析科技有限责任公司利用 LA-ICP-MS 同时分析完成。有关详细的仪器参数和分析流程见 [Zong et al. \(2017\)](#)。GeolasPro 激光剥蚀系统由 COMPexPro 102 ArF 193 nm 准分子激光器和 MicroLas 光学系统组成, ICP-MS 型号为 Agilent 7900。激光剥蚀过程中, 采用氦气作载气, 氩气为补偿气以调节灵敏度。二者在进入 ICP 之前通过一个 T 型接头混合, 激光剥蚀系统配置有信号平滑装置([Hu et al., 2015](#))。本次分析的激光束斑和频率分别为 $\times\times\mu\text{m}$ 和 $\times\times\text{Hz}$ 。U-Pb 同位素定年和微量元素含量处理中, 使用锆石标准 91500 和玻璃标准物质 NIST610 作为外标, 分别进行同位素和微量元素分馏校正([Hu et al., 2021](#))。每个时间分辨分析数据包括 20-30 秒空白信号和 50 秒样品信号。对分析数据的离线处理(包括对样品和空白信号的选择、仪器灵敏度漂移校正、元素含量及 U-Pb 同位素比值和年龄计算)采用软件 ICPMSDataCal ([Liu et al., 2008; Liu et al., 2010](#))完成。对锆石样品的 U-Pb 年龄谱和图绘制和年龄加权平均计算采用 Isoplot/Ex_ver3 ([Ludwig, 2003](#)) 完成。

注: 岩石薄片锆石 U-Pb 同位素定年参考 [Zong et al. \(2010\)](#)。 $\times\times$ 代表实验中激光束斑和频率的数值, 根据具体实验条件填写。锆石常规分析条件为 $32\mu\text{m}$ 和 5Hz 。

1.2 In-situ U-Pb dating and trace element analysis of zircon by LA-ICP-MS (91500)

U-Pb dating and trace element analysis of zircon were simultaneously conducted by LA-ICP-MS at the Wuhan SampleSolution Analytical Technology Co., Ltd., Wuhan, China. Detailed operating conditions for the laser ablation system, the ICP-MS instrument, and data reduction are the same as described by [Zong et al. \(2017\)](#). Laser sampling was performed using a GeolasPro laser ablation system that consists of a COMPexPro 102 ArF excimer laser (wavelength of 193 nm and maximum energy of 200 mJ) and a MicroLas optical system. An Agilent 7900 ICP-MS instrument was used to acquire ion-signal intensities. Helium was applied as a carrier gas. Argon was used as the make-up gas and mixed with the carrier gas via a T-connector before entering the ICP. A “wire” signal smoothing device is included in this laser ablation system ([Hu et al., 2015](#)). The spot size and frequency of the laser were set to $\times\times\mu\text{m}$ and $\times\times\text{Hz}$, respectively, in this study. Zircon 91500 and glass NIST610 served as external standards for U-Pb dating and trace element calibration ([Hu et al., 2021](#)), respectively. Each analysis incorporated a background acquisition of approximately 20–30 seconds, followed by 50 seconds of data acquisition from the sample. An Excel-based software, ICPMSDataCal, was used to perform off-line selection and integration of background and analyzed signals, time-drift correction, and quantitative calibration for trace element analysis and U-Pb dating



(Liu et al., 2008; Liu et al., 2010). Concordia diagrams and weighted mean calculations were conducted using Isoplot/Ex_ver3 (Ludwig, 2003).

Note: The real in-situ U-Pb dating of zircon in thin-section can cite reference of Zong et al. (2010).

References

- Zong, K.Q., Klemd, R., Yuan, Y., He, Z.Y., Guo, J.L., Shi, X.L., Liu, Y.S., Hu, Z.C., Zhang, Z.M., 2017. The assembly of Rodinia: The correlation of early Neoproterozoic (ca. 900 Ma) high-grade metamorphism and continental arc formation in the southern Beishan Orogen, southern Central Asian Orogenic Belt (CAOB). *Precambrian Research*, 290, 32–48.
- Hu, Z.C., Zhang, W., Liu, Y.S., Gao, S., Li, M., Zong, K.Q., Chen, H.H., Hu, S.H., 2015. “Wave” signal-smoothing and mercury-removing device for laser ablation quadrupole and multiple collector ICPMS analysis: application to lead isotope analysis. *Analytical Chemistry*, 87(2), 1152–1157.
- Liu, Y.S., Hu, Z.C., Gao, S., Günther, D., Xu, J., Gao, C.G. and Chen, H.H., 2008. In situ analysis of major and trace elements of anhydrous minerals by LA-ICP-MS without applying an internal standard. *Chemical Geology*, 257(1-2): 34-43.
- Liu, Y.S., Gao, S., Hu, Z.C., Gao, C.G., Zong, K.Q. and Wang, D.B., 2010. Continental and oceanic crust recycling-induced melt-peridotite interactions in the Trans-North China Orogen: U-Pb dating, Hf isotopes and trace elements in zircons of mantle xenoliths. *Journal of Petrology*, 51(1–2): 537–571.
- Ludwig, K.R., 2003. ISOPLOT 3.00: A Geochronological Toolkit for Microsoft Excel. Berkeley Geochronology Center, California, Berkeley, 39 pp.
- Zong, K.Q., Liu, Y.S., Gao, C.G., Hu, Z.C., Gao, S., Gong, H.J., 2010. In situ U-Pb dating and trace element analysis of zircons in thin sections of eclogite: Refining constraints on the ultra-high pressure metamorphism of the Sulu terrane, China. *Chem. Geol.* 269: 237-251.