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5.6 Scheme for Pb isotope ratio analyses using MC-ICP-MS

All chemical preparations were performed on class 100 work benches within a class 1000 overpressured clean laboratory. **Sample digestion:** (1) Sample powder (200 mesh) was placed in an oven at 105 °C for drying of 12 hours; (2) 50–200 mg sample powder was accurately weighed and placed in a Teflon bomb; (3) 1-3 ml HNO₃ and 1-3 ml HF were added into the Teflon bomb; (4) Teflon bomb was putted in a stainless steel pressure jacket and heated to 190 °C in an oven for >24 hours; (5) After cooling, the Teflon bomb was opened and placed on a hotplate at 140 °C and evaporated to incipient dryness, and then 1 ml HNO₃ was added and evaporated to dryness again; (6) The sample was dissolved in 1.0 mL of 1.0 M HBr. **Column chemistry:** After centrifugation, the supernatant solution was loaded into an ion-exchange column packed with AG resin. After complete draining of the sample solution, columns were rinsed with 1.0 M HBr to remove undesirable matrix elements. Finally, the Pb fraction was eluted using 6.0 M HCl and gently evaporated to dryness prior to mass spectrometric measurement.

Pb isotope analyses were performed on a Neptune Plus MC-ICP-MS (Thermo Fisher Scientific, Dreieich, Germany) at the Wuhan Sample Solution Analytical Technology Co., Ltd, Hubei, China. The Neptune Plus, a double focusing MC-ICP-MS, was equipped with seven fixed electron multiplier ICs and nine Faraday cups fitted with 10^{11} Ω resistors. The faraday collector configuration of the mass system was composed of an array to monitor $^{204}(\text{Pb}+\text{Hg})$ 、 ^{206}Pb 、 ^{207}Pb 、 ^{208}Pb 、 ^{203}Tl 、 ^{205}Tl and ^{202}Hg . The large dry interface pump (120 m³ hr⁻¹ pumping speed) and newly designed X skimmer cone and Jet sample cone were used to increase the instrumental sensitivity. A Pb single element solution from Alfa (Alfa Aesar, Karlsruhe, Germany) was used to optimize instrument operating parameters. An aliquot of the international standard solution of 200 μg L⁻¹ NIST 981 was

regularly used for evaluating the reproducibility and accuracy of the instrument. Typically, the signal intensities of $^{208}\text{Pb}^+$ in NIST 981 were $> \sim 6.0$ V. The Pb isotopic data were acquired in the static mode at low resolution. The routine data acquisition consisted of ten blocks of 10 cycles (4.194 s of integration time per cycle). The total time of one measurement lasted about 7 minutes.

The exponential law, which was initially developed for TIMS measurement (Russell et al. 1978) and remains the most widely accepted and utilized with MC-ICP-MS, was used to assess the instrumental mass discrimination in this study. Mass discrimination correction was carried out via normalization to a $^{205}\text{Tl}/^{203}\text{Tl}$ ratio of 2.38714 (the certified value of NIST SRM 997). All data reduction for the MC-ICP-MS analysis of Pb isotope ratios was conducted using “Iso-Compass” software (Zhang et al. 2020). Because of the difference in mass bias behaviors between Pb and Tl, all measured $^{20x}\text{Pb}/^{204}\text{Pb}$ ratios of unknown samples were normalized to the well-accepted NIST 981 values of $^{208}\text{Pb}/^{204}\text{Pb}=36.7262\pm 31$, $^{207}\text{Pb}/^{204}\text{Pb}=15.5000\pm 13$, $^{206}\text{Pb}/^{204}\text{Pb}=16.9416\pm 13$ ($n=119$, Baker et al. 2004). One NIST 981 standard was measured every seven samples analyzed. Analyses of NIST 981 standard yielded external precisions of 0.03% (2RSD) for $^{20x}\text{Pb}/^{204}\text{Pb}$ ratios. In addition, the USGS reference materials BCR-2 (basalt), (maybe different, the test report shall prevail) yielded results of $^{208}\text{Pb}/^{204}\text{Pb}=38.736\pm 17$, $^{207}\text{Pb}/^{204}\text{Pb}=15.628\pm 3$, $^{206}\text{Pb}/^{204}\text{Pb}=18.756\pm 10$ (2SD, $n=22$), respectively, which are identical within an error of 0.03% to their published values (Zhang and Hu 2020).

References

- Russell, W.A., Papanastassiou, D.A., Tombrello, T.A., (1978). Ca isotope fractionation on the earth and other solar system materials. *Geochim. Cosmochim. Acta*, 42 (8), 1075–1090.
- Zhang W., Hu Z.C. (2020). Estimation of isotopic reference values for pure materials and geological reference materials. *At. Spectrosc.* 2020, 41 (3), 93–102.
- Zhang W., Hu Z.C., Liu Y.S. (2020). Iso-Compass: new freeware software for isotopic data reduction of LA-MC-ICP-MS. *J. Anal. At. Spectrom.*, 2020, 35, 1087–1096.
- Baker, J., Peate, D., Waight, T., Meyzen, C. (2004). Pb isotopic analysis of standards and samples using a ^{207}Pb – ^{204}Pb double spike and thallium to correct for mass bias with a double-focusing mc-icp-ms. *Chemical Geology*, 211(3–4), 275–303