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1.1 In-situ U-Pb dating and trace element analysis of zircon by LA-ICP-MS (Tanz)

U-Pb dating and trace element analysis of zircon were simultaneously conducted by LA-ICP-MS at the Wuhan SampleSolution Analytical Technology Co., Ltd., Wuhan, China. Detailed operating conditions for the laser ablation system, the ICP-MS instrument, and data reduction are the same as described by [Zong et al. \(2017\)](#). Laser sampling was performed using a GeolasPro laser ablation system that consists of a COMPexPro 102 ArF excimer laser (wavelength of 193 nm and maximum energy of 200 mJ) and a MicroLas optical system. An Agilent 7900 ICP-MS instrument was used to acquire ion-signal intensities. Helium was applied as a carrier gas. Argon was used as the make-up gas and mixed with the carrier gas via a T-connector before entering the ICP. A “wire” signal smoothing device is included in this laser ablation system ([Hu et al., 2015](#)). The spot size and frequency of the laser were set to $\times\times\mu\text{m}$ and $\times\times\text{Hz}$, respectively, in this study. Zircon Tanz and glass NIST610 served as external standards for U-Pb dating and trace element calibration ([Hu et al., 2021](#)),

respectively. Each analysis incorporated a background acquisition of approximately 20–30 seconds, followed by 50 seconds of data acquisition from the sample. An Excel-based software, ICPMSDataCal, was used to perform off-line selection and integration of background and analyzed signals, time-drift correction, and quantitative calibration for trace element analysis and U-Pb dating (Liu et al., 2008; Liu et al., 2010). Concordia diagrams and weighted mean calculations were conducted using Isoplot/Ex_ver3 (Ludwig, 2003).

Note: The real in-situ U-Pb dating of zircon in thin-section can cite reference of Zong et al. (2010).

References

- Hu, Z.C., Li, X.H., Luo, T., Zhang, W., Crowley, J., Li, Q.L., Ling, X.X., Yang, C., Li, Y., Feng, L.P., Xia, X.P., Zhang, S.B., Wang, Z.C., Guo, J.L., Xu, L., Lin, J., Liu, X.M., Bao, Z.A., Liu, Y.S., Zong, K.Q., Chen, W., Hu, S. H., 2021. Tanz zircon megacrysts: a new zircon reference material for the microbeam determination of U–Pb ages and Zr–O isotopes. *Journal of Analytical Atomic Spectrometry*, 36(12), 2715–2734.
- Hu, Z.C., Zhang, W., Liu, Y.S., Gao, S., Li, M., Zong, K.Q., Chen, H.H., Hu, S.H., 2015. “Wave” signal-smoothing and mercury-removing device for laser ablation quadrupole and multiple collector ICPMS analysis: application to lead isotope analysis. *Analytical Chemistry*, 87(2), 1152–1157.
- Liu, Y.S., Gao, S., Hu, Z.C., Gao, C.G., Zong, K.Q. and Wang, D.B., 2010. Continental and oceanic crust recycling-induced melt-peridotite interactions in the Trans-North China Orogen: U-Pb dating, Hf isotopes and trace elements in zircons of mantle xenoliths. *Journal of Petrology*, 51(1–2): 537–571.
- Liu, Y.S., Hu, Z.C., Gao, S., Günther, D., Xu, J., Gao, C.G. and Chen, H.H., 2008. In situ analysis of major and trace elements of anhydrous minerals by LA-ICP-MS without applying an internal standard. *Chemical Geology*, 257(1-2): 34–43.
- Ludwig, K.R., 2003. ISOPLOT 3.00: A Geochronological Toolkit for Microsoft Excel. Berkeley Geochronology Center, California, Berkeley, 39 pp.
- Zong, K.Q., Klemd, R., Yuan, Y., He, Z.Y., Guo, J.L., Shi, X.L., Liu, Y.S., Hu, Z.C., Zhang, Z.M., 2017. The assembly of Rodinia: The correlation of early Neoproterozoic (ca. 900 Ma) high-grade metamorphism and continental arc formation in the southern Beishan Orogen, southern Central Asian Orogenic Belt (CAOB). *Precambrian Research*, 290, 32–48.
- Zong, K.Q., Liu, Y.S., Gao, C.G., Hu, Z.C., Gao, S., Gong, H.J., 2010. In situ U-Pb dating and trace element analysis of zircons in thin sections of eclogite: Refining constraints on the ultra-high pressure metamorphism of the Sulu terrane, China. *Chem. Geol.* 269: 237–251.