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### 5.12 Mg isotope analysis using MC-ICP-MS

All chemical preparations were performed on class 100 work benches within a class 1000 overpressured clean laboratory. **Sample digestion :** Based on the MgO content of the sample, a portion of each sample powder containing 20  $\mu\text{g}$  of Mg was weighed into Teflon bombs. Then, 2ml  $\text{HNO}_3$ -HF mixed acid (1:1, v/v) was added, and the sealed bomb was placed in an electric oven at 190°C for 48 hours to digest the sample. After cooling, the bomb was opened and placed on a hotplate to evaporated to dryness. This was followed by adding 1 ml  $\text{HNO}_3$  and evaporating to dryness. This step was repeated twice to completely remove fluoride. After that, the residue was re-dissolved by adding 1 ml of  $\text{HNO}_3$  and 2 ml of ultra-pure water. The bomb was resealed and placed in an electric oven at 190°C for 12 hours to obtain a clear digestion solution. The final solution was evaporated and re-dissolved in 200  $\mu\text{l}$  1.5N  $\text{HNO}_3$  for chemical separation.

**Magnesium separation:**The chemical separation procedure for geological samples followed method outlined by Wang et al. (2022), and involved three sequential column separation steps. Step one: The digested sample was loaded into a column filled with AG50W-X8 resin (200~400 mesh), and subsequently rinsed with 20 ml of 1.5N  $\text{HNO}_3$ . All fractions from both loading and rinsing were collected for the quantitative recovery of Mg. The residual matrix elements (such as Al, Fe, Ca) are then eluted from the resin using 6N HCl. Step two: The Mg solution obtained in step one was evaporated and loaded to the same column which was pre-cleared. Matrix elements are washed with 1N  $\text{HNO}_3$ , and the Mg fraction was collected through the elution with 17 ml 1N  $\text{HNO}_3$ . Step three: Repeat step two to further improve the purity of the Mg solution.

After the three-column chemical purification process, all the interference elements were effectively separated from Mg. The recovery of Mg obtained was 100%±3%, and the total blank of the entire process was less than 12 ng.

**Mass spectrometry measurements:**Mg isotopic measurement was performed on a Thermo Fisher Scientific Neptune Plus MC-ICP-MS at Wuhan Sample Solution Analytical Technology Co., Ltd. Measurements were performed in low resolution and sample solution were introduced into the plasma at 200pb Mg in 2%  $\text{HNO}_3$  (m/m). The sample-standard bracketing (SSB) method was used to correct the instrumental mass bias during Mg isotopic analysis. All Mg isotope data are measured

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### Wuhan SampleSolution Analytical Technology Co., Ltd

relative to the GSB Mg standard, and then converted to DSM3 based results according to the relative value between GSB Mg and DSM3 proposed by Bao et al., 2019. The converted formula is as follows:

$$d^i\text{Mg} (\text{‰}) = [(^i\text{Mg}/^{24}\text{Mg})_{\text{sample}} / (^i\text{Mg}/^{24}\text{Mg})_{\text{GSB Mg}} - 1] \times 10^3 \quad (\text{where } i = 25 \text{ or } 26)$$

$$\delta^{26/24}\text{Mg}_{\text{DSM3}} = \delta^{26/24}\text{Mg}_{\text{smp-GSB Mg}} - 2.05$$

$$\delta^{25/24}\text{Mg}_{\text{DSM3}} = \delta^{25/24}\text{Mg}_{\text{smp-GSB Mg}} - 1.06$$

Each batch of samples includes at least two international geological reference standards to monitor the performance of instrumental measurement and chemical separation, ensuring the accuracy and precision of the results. Repeated measurements of the geological standards (e.g., BCR-2, AGV-2, BHVO-2 and GSP-2, the test report shall prevail) in our laboratory were consistent with previously reported data, and the long-term reproducibility of  $\delta^{26}\text{Mg}$  was better than 0.06‰.

#### References

Bao, Z.-A., K.-J. Huang, T.-Z. Huang, B. Shen, C.-L. Zong, K.-Y. Chen, and H.-L. Yuan, 2019, Precise magnesium isotope analyses of high-K and low-Mg rocks by MC-ICP-MS: *Journal of Analytical Atomic Spectrometry*, v. 34, p. 940-953.

Wang, Y., Y. Zhang, X. Li, S. Ke, A. Sun, R. Yang, W. Yang, and Y. He, 2022, Purification of Mg from extremely low-Mg felsic rocks for isotopic ratio determination by MC-ICP-MS: *Journal of Analytical Atomic Spectrometry*, v. 37, p. 497-507.