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### 1.9 In-situ U-Pb dating and trace element analysis of Apatite by LA-ICP-MS

U-Pb dating and trace element analysis of apatite were simultaneously conducted by LA-ICP-MS at the Wuhan SampleSolution Analytical Technology Co., Ltd., Wuhan, China. Laser sampling was performed using a GeolasPro laser ablation system that consists of a COMPexPro 102 ArF excimer laser (wavelength of 193 nm and maximum energy of 200 mJ) and a MicroLas optical system. An Agilent 7900 ICP-MS instrument was used to acquire ion-signal intensities. Helium was applied as a carrier gas. Argon was used as the make-up gas and mixed with the carrier gas via a T-connector before entering the ICP. A “wire” signal smoothing device is included in this laser ablation system (Hu et al., 2015). The spot size and frequency of the laser were set to  $\times\times\mu\text{m}$  and  $\times\times\text{Hz}$ , respectively. Apatite MAD (Thomson et al., 2012) was used as an external standard to correct the Pb/U fractionation and instrumental mass discrimination, The trace element compositions of apatite were calibrated against NIST 612 glass as an external calibration (Liu et al., 2008). Apatite Sumé was analyzed as an unknown. The obtained weighted average  $^{206}\text{Pb}/^{238}\text{U}$  ages in  $\times$  analyses of Sumé are  $***\text{Ma}$  in this study, which are consistent with the reference age of  $568.0\pm 3.1\text{Ma}$  (Lana et al. 2022). Each analysis incorporated a background acquisition of approximately 20-30 s followed by 50 s of data acquisition from the sample. A software Iolite 4.0 was used to perform off-line selection and

integration of background and analyzed signals, time-drift correction, and quantitative calibration for trace element analysis and U-Pb dating (Paton et al. 2011). Tera-wasserburg diagrams of U-Pb ages and  $^{207}\text{Pb}$  corrected age weighted average calculations were made using Isoplot/Ex\_ver3 (Ludwig, 2003).

### References

- Hu, Z.C., Zhang, W., Liu, Y.S., Gao, S., Li, M., Zong, K.Q., Chen, H.H., Hu, S.H., 2015. “Wave” signal-smoothing and mercury-removing device for laser ablation quadrupole and multiple collector ICPMS analysis: application to lead isotope analysis. *Analytical Chemistry*, 87(2), 1152–1157.
- Thomson, S. N., G. E. Gehrels, J. Ruiz, and R. Buchwaldt (2012), Routine low-damage apatite U-Pb dating using laserablation–multicollector–ICPMS, *Geochem. Geophys. Geosyst.*, 13, Q0AA21
- Liu, Y.S., Hu, Z.C., Gao, S., Günther, D., Xu, J., Gao, C.G. and Chen, H.H., 2008. In situ analysis of major and trace elements of anhydrous minerals by LA-ICP-MS without applying an internal standard. *Chemical Geology*, 257(1-2): 34-43.
- Lana, C., Gonçalves, G.O., Mazoz, A., Buick, I., Kamo, S., Scholz, R., Wang, H., Moreira, H., Babinski, M. and Queiroga, G. (2022), Assessing the U-Pb, Sm-Nd and Sr-Sr Isotopic Compositions of the Sumé Apatite as a Reference Material for LA-ICP-MS Analysis. *Geostand Geoanal Res*, 46: 71-95.
- Paton, C., Hellstrom, J., Paul, B., Woodhead, J. and Hergt, J. (2011) Iolite: Freeware for the visualisation and processing of mass spectrometric data. *Journal of Analytical Atomic Spectrometry*. 26, 2508–2518 .
- Ludwig, K.R., 2003. ISOPLOT 3.00: A Geochronological Toolkit for Microsoft Excel. Berkeley Geochronology Center, California, Berkeley, 39 pp.