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1.7 In-situ U-Pb dating and trace element analysis of xenotime by LA-ICP-MS

U-Pb dating and trace element analyses of xenotime were conducted at the Wuhan SampleSolution Analytical Technology Co., Ltd., Wuhan, China. Laser sampling was performed using a GeolasPro HD laser ablation system that consists of a COMPexPro 102 ArF excimer laser (wavelength of 193 nm and maximum energy of 200 mJ) and a MicroLas optical system. An Agilent 7900 ICP-MS instrument was used to acquire ion-signal intensities. All analyses were performed with a laser spot size of $\times\mu\text{m}$, a repetition rate of $\times\text{ Hz}$ and a fluence of 5 J/cm^2 in this study. The Helium was used as the carrier gas in the ablation cell and merged with argon (makeup gas) behind the ablation cell (Günther and A. Heinrich, 1999; Luo et al., 2018a). A signal-smoothing and mercury-removing device was used in this laser ablation system to obtain smooth signals and reduce the mercury signal (Hu et al., 2014). A small amount of (4.1 mg min^{-1}) water vapor was added before the ablation cell to improve the analytical accuracy and precision (Luo et al., 2018b). Each single-spot analysis consisted of 20 seconds of background signal acquisition followed by 50 seconds of ablation. Zircon 91500 (Wiedenbeck et al., 1995) was used as an external standard to correct the

Pb/U fractionation and instrumental mass discrimination, and xenotime BS-1 was analyzed as an unknown. The obtained weighted average $^{206}\text{Pb}/^{238}\text{U}$ ages in x analyses of BS-1 are ***Ma in this study, which are consistent with the reference age of 508.9 ± 0.3 Ma (Fletcher et al., 2004). The trace element compositions of xenotime were calibrated against NIST 610 glass as an external calibration (Liu et al., 2010a). Off-line selection and integration of background and analyte signals, time-drift correction, and quantitative calibration for trace element analyses and U-Pb dating were performed by ICPMSDataCal (Liu et al., 2010b). Concordia diagrams and weighted mean calculations were made using Isoplot/Ex_ver3 (Ludwig, 2003).

References

- Günther, D., A. Heinrich, C., 1999. Enhanced sensitivity in laser ablation-ICP mass spectrometry using helium-argon mixtures as aerosol carrier. *Journal of Analytical Atomic Spectrometry*, 14(9): 1363-1368.
- Luo, T. et al., 2018a. Reassessment of the influence of carrier gases He and Ar on signal intensities in 193nm excimer LA-ICP-MS analysis. *Journal of Analytical Atomic Spectrometry*, 33(10): 1655-1663.
- Hu, Z.C., Zhang, W., Liu, Y.S., Gao, S., Li, M., Zong, K.Q., Chen, H.H., Hu, S.H., 2015. "Wave" signal-smoothing and mercury-removing device for laser ablation quadrupole and multiple collector ICPMS analysis: application to lead isotope analysis. *Analytical Chemistry*, 87(2), 1152-1157.
- Luo, T. et al., 2018b. Water Vapor-Assisted "Universal" Nonmatrix-Matched Analytical Method for the in Situ U-Pb Dating of Zircon, Monazite, Titanite, and Xenotime by Laser Ablation-Inductively Coupled Plasma Mass Spectrometry. *Analytical chemistry*, 90(15): 9016-9024.
- Wiedenbeck, M. et al., 1995. Three natural zircon standards for U-Th-Pb, Lu-Hf, trace element and REE analysis. *Geostandards Newsletter*, 19(1): 1-23.
- Fletcher, I.R., McNaughton, N.J., Aleinikoff, J.A., Rasmussen, B., Kamo, S.L., 2004. Improved calibration procedures and new standards for U-Pb and Th-Pb dating of Phanerozoic xenotime by ion microprobe. *Chemical Geology*, 209(3-4): 295-314.
- Liu, Y.S., Hu, Z.C., Zong, K.Q., Gao, C.G., Gao, S., Xu, J. and Chen, H.H., 2010a. Reappraisal and refinement of zircon U-Pb isotope and trace element analyses by LA-ICP-MS. *Chinese Science Bulletin*, 55(15): 1535-1546.
- Liu, Y.S., Gao, S., Hu, Z.C., Gao, C.G., Zong, K.Q. and Wang, D.B., 2010b. Continental and oceanic crust recycling-induced melt-peridotite interactions in the Trans-North China Orogen: U-Pb dating, Hf isotopes and trace elements in zircons of mantle xenoliths. *Journal of Petrology*, 51(1&2): 537-571.
- Ludwig, K.R., 2003. ISOPLOT 3.00: A Geochronological Toolkit for Microsoft Excel. Berkeley Geochronology Center, California, Berkeley, 39 pp.