LA-ICP-MS/LA-MC-ICP-MS 数据处理软件

# **ICPMSDataCal**

**(V11.4)**

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刘勇胜 (yshliu@hotmail.com)

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1. Liu Y.S., Hu Z.C., Gao S., Günther D., Xu J., Gao C.G. and Chen H.H., 2008. In situ analysis of major and trace elements of anhydrous minerals by LA-ICP-MS without applying an internal standard. *Chemical Geology*, 257(1-2): 34-43.
2. Chen L., Liu Y.S., Hu Z.C., Gao S., Zong K.Q. and Chen H.H., 2011. Accurate determinations of fifty-four major and trace elements in carbonate by LA-ICP-MS using normalization strategy of bulk components as 100%. *Chemical Geology*, 284(3-4): 283-295.

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1. Liu Y.S., Gao S., Hu Z.C., Gao C.G., Zong K.Q. and Wang D.B., 2010. Continental and oceanic crust recycling-induced melt-peridotite interactions in the Trans-North China Orogen: U-Pb dating, Hf isotopes and trace elements in zircons from mantle xenoliths. *Journal of Petrology*, 51(1&2): 537-571.
2. Liu Y.S., Hu Z.C., Zong K.Q., Gao C.G., Gao S., Xu J. and Chen H.H., 2010b. Reappraisal and refinement of zircon U-Pb isotope and trace element analyses by LA-ICP-MS. *Chinese Science Bulletin*, 55(15): 1535-1546.
3. Hu Z.C., Liu Y.S., Gao S., Liu W.G., Yang L., Zhang W., Tong X.R., Lin L., Zong K.Q., Li M., Chen H.H. and Zhou L., 2012. Improved in situ Hf isotope ratio analysis of zircon using newly designed X skimmer cone and Jet sample cone in combination with the addition of nitrogen by laser ablation multiple collector ICP-MS. *J. Anal. At. Spectrom.*, 27: 1391-1399

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1. Jie Lin, **Yongsheng Liu**, Yueheng Yang, Zhaochu Hu, 2016. Calibration and correction of LA-ICP-MS and LA-MC-ICP-MS analyses for element contents and isotopic ratios. ***Solid Earth Sciences***, 1(1): 5–27.
2. \*Jie Lin, **Yongsheng Liu**, Zhaochu Hu, Lu Yang, Kang Chen, Haihong Chen, Keqing Zong and Shan Gao, 2016. Accurate determination of lithium isotope ratios by MC-ICP-MS without strict matrix-matching by using a novel washing method. ***J. Anal. At. Spectrom***, **31**, 390 – 397.
3. \*Jie Lin, Yongsheng Liu, Haihong Chen, Lian Zhou, Zhaochu Hu, Shan Gao. Review of High-Precision Sr Isotope Analyses of Low-Sr Geological Samples. ***Journal of Earth Science***, 26(5): 763-774.

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This screenshot shows the 'IR 同位素比值校正方法及设置' (IR Isotope Ratio Correction Method and Settings) dialog box for Sr isotope analysis. The '分馏因子平均值 (FF)' (Fractionation Factor Average) section is circled in red, showing values for '88Sr/86Sr' (-1.3205) and '87Rb/85Sr' (-1.4678). The '同位素标准列表' (Isotope Standard List) section shows 'JUN12C01' selected. The '分馏和漂移校正' (Fractionation and Drift Correction) section has checkboxes for '采用外标校正' (Use external standard correction) and '时间漂移校正' (Time drift correction). The '分馏因子计算方法 (Sr)' (Fractionation Factor Calculation Method (Sr)) section has radio buttons for '采用实时的分馏因子' (Use real-time fractionation factor), '采用分馏因子平均值' (Use fractionation factor average), and '采用外标的分馏因子' (Use external standard fractionation factor). The '分馏因子计算方法 (Rb | StHs6/80-G)' (Fractionation Factor Calculation Method (Rb | StHs6/80-G)) section has radio buttons for 'Rb标准独立测定' (Rb standard independent determination), '利用外标计算值' (Use external standard calculation value), and 'FF' (Fractionation Factor). The 'Ignore 87Rb interferences on 87Sr' checkbox is unchecked.

This screenshot shows the 'IR 同位素比值校正方法及设置' (IR Isotope Ratio Correction Method and Settings) dialog box for Nd isotope analysis. The '分馏因子平均值 (FF)' (Fractionation Factor Average) section is circled in red, showing values for '146Nd/144Nd' (-1.4291) and '147Sm/149Sm' (-1.2388). The '同位素标准列表' (Isotope Standard List) section shows 'DEC06A01' selected. The '分馏和漂移校正' (Fractionation and Drift Correction) section has checkboxes for '采用外标校正' (Use external standard correction) and '时间漂移校正' (Time drift correction). The '分馏因子计算方法 (Nd)' (Fractionation Factor Calculation Method (Nd)) section has radio buttons for '采用实时的分馏因子' (Use real-time fractionation factor), '采用分馏因子平均值' (Use fractionation factor average), and '采用外标的分馏因子' (Use external standard fractionation factor). The '分馏因子计算方法 (Sm)' (Fractionation Factor Calculation Method (Sm)) section has radio buttons for '采用样品实测值' (Use sample measured value) and 'Sm标准独立测定' (Sm standard independent determination). The 'Ignore 144Sm interferences on 144Nd' checkbox is unchecked.

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Chen, L., Liu, Y., Hu, Z., Gao, S., Zong, K. and Chen, H., 2011. Accurate determinations of fifty-four major and trace elements in carbonate by LA-ICP-MS using normalization strategy of bulk components as 100%. *Chemical Geology*, 284(3-4): 283-295.

Liu, Y.S., Hu, Z.C., Gao, S., Günther, D., Xu, J., Gao, C.G. & Chen, H.H., 2008. In situ analysis of major and trace elements of anhydrous minerals by LA-ICP-MS without applying an internal standard. *Chemical Geology* **257**, 34-43.

Liu, Y.S., Gao, S., Hu, Z.C., Gao, C.G., Zong, K.Q. and Wang, D.B., 2009. Continental and oceanic crust recycling-induced melt-peridotite interactions in the Trans-North China Orogen: U-Pb dating, Hf isotopes and trace elements in zircons from mantle xenoliths. *Journal of Petrology*, 51(1&2): 537-571.

Liu, Y.S., Hu, Z.C., Zong, K.Q., Gao, C.G., Gao, S., Xu, J., and Chen, H.H., 2010. Reappraisal and refinement of zircon U-Pb isotope and trace element analyses by LA-ICP-MS. *Chinese Sci. Bull.*, 55(15): 1535-1546.

Heinrich, C. A., Pettke, T., Halter, W. E., Aigner-Torres, M., Audetat, A., Günther, D., Hattendorf, B., Bleiner, D., Guillong, M. & Horn, I. (2003). Quantitative multi-element analysis of minerals, fluid and melt inclusions by laser-ablation inductively-coupled-plasma mass-spectrometry. *Geochimica et Cosmochimica Acta* **67**, 3473-3497.

Heinrich, C. A., Ryan, C. G., Mernagh, T. P. & Eadington, P. J. (1992). Segregation of ore metals between magmatic brine and vapor-A fluid inclusion study using pixel microanalysis. *Econ. Geol.* **87**, 1566-1583.

Zajacz, Z. & Halter, W. (2007). LA-ICPMS analyses of silicate melt inclusions in co-precipitated minerals: Quantification, data analysis and mineral/melt partitioning. *Geochimica et Cosmochimica Acta* **71**, 1021-1040.

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## (1) U-Pb dating and trace element analyses of zircon by LA-ICP-MS

U-Pb dating and trace element analyses of zircon were conducted synchronously by LA-ICP-MS at the State Key Laboratory of Geological Processes and Mineral Resources, China University of Geosciences, Wuhan. Detailed operating conditions for the laser ablation system and the ICP-MS instrument and data reduction are the same as description by [Liu \*et al.\* \(2010a; 2008a; 2010b\)](#). Laser sampling was performed using a GeoLas 2005. An Agilent 7500a ICP-MS instrument was used to acquire ion-signal intensities. A “wire” signal smoothing device is included in this laser ablation system, by which smooth signals are produced even at very low laser repetition rates down to 1 Hz ([Hu \*et al.\*, 2012b](#)). Helium was applied as a carrier gas. Argon was used as the make-up gas and mixed with the carrier gas via a T-connector before entering the ICP. Nitrogen was added into the central gas flow (Ar+He) of the Ar plasma to decrease the detection limit and improve precision ([Hu \*et al.\*, 2008a; Liu \*et al.\*, 2010b](#)). Each analysis incorporated a background acquisition of approximately 20-30 s (gas blank) followed by 50 s of data acquisition from the sample. The Agilent Chemstation was utilized for the acquisition of each individual analysis. Off-line selection and integration of background and analyte signals, and time-drift correction and

quantitative calibration for trace element analyses and U-Pb dating were performed by *ICPMSDataCal* (Liu et al., 2010a; Liu et al., 2008a).

Zircon 91500 was used as external standard for U-Pb dating, and was analyzed twice every 5 analyses. Time-dependent drifts of U-Th-Pb isotopic ratios were corrected using a linear interpolation (with time) for every five analyses according to the variations of 91500 (i.e., 2 zircon 91500 + 5 samples + 2 zircon 91500)(Liu et al., 2010a). Preferred U-Th-Pb isotopic ratios used for 91500 are from Wiedenbeck et al. (1995). Uncertainty of preferred values for the external standard 91500 was propagated to the ultimate results of the samples. Concordia diagrams and weighted mean calculations were made using Isoplot/Ex\_ver3 (Ludwig, 2003). Trace element compositions of zircons were calibrated against multiple-reference materials (BCR-2G and BIR-1G) combined with internal standardization (Liu et al., 2010a) (或者 without applying internal standardization (Liu et al., 2008a), 具体分析采用的方法会因个人处理数据的方式不同而有变化, 如不清楚可咨询实验室管理员) . The preferred values of element concentrations for the USGS reference glasses are from the GeoReM database (<http://georem.mpch-mainz.gwdg.de/>).

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\*Zong KQ, Liu YS, Gao CG, Hu ZC, Gao S, Gong HJ, 2010. In situ U-Pb

dating and trace element analysis of zircons in thin sections of eclogite: Refining constraints on the ultra-high pressure metamorphism of the Sulu terrane, China. *Chem. Geol.* 269: 237-251. DOI: 10.1016/j.chemgeo.2009.09.021.

## (2) Trace element analyses by LA-ICP-MS

Major and trace element analyses were conducted by LA-ICP-MS at the State Key Laboratory of Geological Processes and Mineral Resources, China University of Geosciences, Wuhan. Detailed operating conditions for the laser ablation system and the ICP-MS instrument and data reduction are the same as description by [Liu \*et al.\* \(2008a\)](#). Laser sampling was performed using a GeoLas 2005. An Agilent 7500a ICP-MS instrument was used to acquire ion-signal intensities. A “wire” signal smoothing device is included in this laser ablation system, by which smooth signals are produced even at very low laser repetition rates down to 1 Hz ([Hu \*et al.\*, 2012b](#)). Helium was applied as a carrier gas. Argon was used as the make-up gas and mixed with the carrier gas via a T-connector before entering the ICP. Nitrogen was added into the central gas flow (Ar+He) of the Ar plasma to decrease the detection limit and improve precision ([Hu \*et al.\*, 2008a](#)). Each analysis incorporated a background acquisition of approximately

20-30 s (gas blank) followed by 50 s of data acquisition from the sample. The Agilent Chemstation was utilized for the acquisition of each individual analysis. Element contents were calibrated against multiple-reference materials (BCR-2G, BIR-1G and BHVO-2G) without applying internal standardization (Liu et al., 2008a). The preferred values of element concentrations for the USGS reference glasses are from the GeoReM database (<http://georem.mpch-mainz.gwdg.de/>). Off-line selection and integration of background and analyte signals, and time-drift correction and quantitative calibration were performed by *ICPMSDataCal* (Liu et al., 2010a; Liu et al., 2008a).

### (3) In situ Hf isotope ratio analysis of zircon by LA-MC-ICP-MS

Experiments were conducted using a Neptune Plus MC-ICP-MS (Thermo Fisher Scientific, Germany) in combination with a Geolas 2005 excimer ArF laser ablation system (Lambda Physik, Göttingen, Germany) that was hosted at the state Key Laboratory of Geological Processes and Mineral Resources, China University of Geosciences in Wuhan. The energy density of laser ablation that was used in this study was  $5.3 \text{ J cm}^{-2}$ . A “wire” signal smoothing device is included in this laser ablation system, by which smooth signals are produced even

at very low laser repetition rates down to 1 Hz (Hu et al., 2012b). Helium was used as the carrier gas within the ablation cell and was merged with argon (makeup gas) after the ablation cell. As demonstrated by our previous study, for the 193 nm laser a consistent 2-fold signal enhancement was achieved in helium than in argon gas (Hu et al., 2008b). We used a simple Y junction downstream from the sample cell to add small amounts of nitrogen ( $4 \text{ ml min}^{-1}$ ) to the argon makeup gas flow (Hu et al., 2008a). Compared to the standard arrangement, the addition of nitrogen in combination with the use of the newly designed X skimmer cone and Jet sample cone in Neptune Plus improved the signal intensity of Hf, Yb and Lu by a factor of 5.3, 4.0 and 2.4, respectively. All data were acquired on zircon in single spot ablation mode at a spot size of  $44 \text{ }\mu\text{m}$  in this study. Each measurement consisted of 20 s of acquisition of the background signal followed by 50 s of ablation signal acquisition. Detailed operating conditions for the laser ablation system and the MC-ICP-MS instrument and analytical method are the same as description by Hu et al. (2012a).

The major limitation to accurate *in situ* zircon Hf isotope determination by LA-MC-ICP-MS is the very large isobaric interference from  $^{176}\text{Yb}$  and, to a much lesser extent  $^{176}\text{Lu}$  on  $^{176}\text{Hf}$  (Woodhead et al.,

2004). The under- or over-estimation of the  $\beta_{Yb}$  value would undoubtedly affect the accurate correction of  $^{176}Yb$  and thus the determined  $^{176}Hf/^{177}Hf$  ratio. We applied the directly obtained  $\beta_{Yb}$  value from the zircon sample itself in real-time (Liu et al., 2010a). The  $^{179}Hf/^{177}Hf$  and  $^{173}Yb/^{171}Yb$  ratios were used to calculate the mass bias of Hf ( $\beta_{Hf}$ ) and Yb ( $\beta_{Yb}$ ), which were normalized to  $^{179}Hf/^{177}Hf = 0.7325$  and  $^{173}Yb/^{171}Yb = 1.1248$  (Blichert-Toft et al., 1997) using an exponential correction for mass bias. Interference of  $^{176}Yb$  on  $^{176}Hf$  was corrected by measuring the interference-free  $^{173}Yb$  isotope and using  $^{176}Yb/^{173}Yb = 0.7876$  (McCulloch et al., 1977) to calculate  $^{176}Yb/^{177}Hf$ . Similarly, the relatively minor interference of  $^{176}Lu$  on  $^{176}Hf$  was corrected by measuring the intensity of the interference-free  $^{175}Lu$  isotope and using the recommended  $^{176}Lu/^{175}Lu = 0.02656$  (Blichert-Toft et al., 1997) to calculate  $^{176}Lu/^{177}Hf$ . We used the mass bias of Yb ( $\beta_{Yb}$ ) to calculate the mass fractionation of Lu because of their similar physicochemical properties. Off-line selection and integration of analyte signals, and mass bias calibrations were performed using ICPMSDataCal (Liu et al., 2010a).

#### (4) Trace element analyses of whole rock

Whole rock samples were crushed in a corundum jaw crusher (to 60

mesh). About 60 g was powdered in an agate ring mill to less than 200 mesh. The samples were then digested by HF + HNO<sub>3</sub> in Teflon bombs and analyzed with an Agilent 7500a ICP-MS at the State Key Laboratory of Geological Processes and Mineral Resources, China University of Geosciences, Wuhan. The detailed sample-digesting procedure for ICP-MS analyses and analytical precision and accuracy for trace elements are the same as description by [Liu et al. \(2008b\)](#).

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\*Zong KQ, Liu YS, Gao CG, Hu ZC, Gao S, Gong HJ, 2010. In situ U-Pb dating and trace element analysis of zircons in thin sections of eclogite: Refining constraints on the ultra-high pressure metamorphism of the Sulu terrane, China. *Chem. Geol.* 269: 237-251.

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## References

- Blichert-Toft J., Chauvel C., Albarède F., 1997. Separation of Hf and Lu for high-precision isotope analysis of rock samples by magnetic sector-multiple collector ICP-MS. *Contributions to Mineralogy and Petrology*, 127(3): 248-260.
- Chen L., Liu Y.S., Hu Z.C., Gao S., Zong K.Q., Chen H.H., 2011. Accurate determinations of fifty-four major and trace elements in carbonate by LA-ICP-MS using normalization strategy of bulk components as 100%. *Chemical Geology*, 284(3-4): 283-295.
- Gao C.G., Liu Y.S., Zong K.Q., Hu Z.C., Gao S., 2010. Microgeochemistry of rutile and zircon in eclogites from the CCSD main hole: Implications for the fluid activity and thermo-history of the UHP metamorphism. *Lithos*, 115(1-4): 51-64.
- Hu Z.C., Gao S., Liu Y.S., Hu S.H., Chen H.H., Yuan H.L., 2008a. Signal enhancement in laser ablation ICP-MS by addition of nitrogen in the central channel gas. *Journal of Analytical Atomic Spectrometry*, 23: 1093-1101.
- Hu Z.C., Liu Y.S., Gao S., Hu S.H., Dietiker R., Günther D., 2008b. A local aerosol extraction strategy for the determination of the aerosol composition in laser ablation inductively coupled plasma mass spectrometry. *Journal of Analytical Atomic Spectrometry*, 23: 1192-1203.
- Hu Z.C., Liu Y.S., Gao S., Liu W.G., Yang L., Zhang W., Tong X.R., Lin L., Zong K.Q., Li M., Chen H.H., Zhou L., 2012a. Improved in situ Hf isotope ratio analysis of zircon using newly designed X skimmer cone and Jet sample cone in combination with the addition of nitrogen by laser ablation multiple collector ICP-MS. *Journal of Analytical Atomic Spectrometry*, 27: 1391-1399.
- Hu Z.C., Liu Y.S., Gao S., Xiao S.Q., Zhao L.S., Günther D., Li M., Zhang W., Zong K.Q., 2012b. A “wire” signal smoothing device for laser ablation inductively coupled plasma mass spectrometry analysis. *Spectrochimica Acta Part B: Atomic Spectroscopy*, 78(0): 50-57.
- Lin J., Liu Y., Tong X., Zhu L., Zhang W., Hu Z., 2017. Improved in situ Li isotopic ratio analysis of silicates by optimizing signal intensity, isotopic ratio stability and intensity matching using ns-LA-MC-ICP-MS. *Journal of Analytical Atomic Spectrometry*: DOI: 10.1039/C6JA00409A.
- Lin J., Liu Y., Yang Y., Hu Z., 2016. Calibration and correction of LA-ICP-MS and LA-MC-ICP-MS analyses for element contents and isotopic ratios. *Solid Earth Sciences*, 1(1): 5-27.
- Liu Y.S., Gao S., Hu Z.C., Gao C.G., Zong K.Q., Wang D.B., 2010a. Continental and oceanic crust recycling-induced melt-peridotite interactions in the Trans-North China Orogen: U-Pb dating,

- Hf isotopes and trace elements in zircons of mantle xenoliths. *Journal of Petrology*, 51(1&2): 537-571.
- Liu Y.S., Hu Z.C., Gao S., Günther D., Xu J., Gao C.G., Chen H.H., 2008a. In situ analysis of major and trace elements of anhydrous minerals by LA-ICP-MS without applying an internal standard. *Chemical Geology*, 257(1-2): 34-43.
- Liu Y.S., Hu Z.C., Yuan H.L., Hu S.H., Cheng H.H., 2007. Volume-optional and low-memory (VOLM) chamber for laser ablation-ICP-MS: application to the analyses of fiber. *Journal of Analytical Atomic Spectrometry*, 22: 582-585.
- Liu Y.S., Hu Z.C., Zong K.Q., Gao C.G., Gao S., Xu J., Chen H.H., 2010b. Reappraisal and refinement of zircon U-Pb isotope and trace element analyses by LA-ICP-MS. *Chinese Science Bulletin*, 55(15): 1535-1546.
- Liu Y.S., Zong K.Q., Kelemen P.B., Gao S., 2008b. Geochemistry and magmatic history of eclogites and ultramafic rocks from the Chinese continental scientific drill hole: Subduction and ultrahigh-pressure metamorphism of lower crustal cumulates. *Chemical Geology*, 247(1-2): 133-153.
- Ludwig K.R., 2000. Decay constant errors in U-Pb concordia-intercept ages. *Chemical Geology*, 166(3-4): 315-318.
- Ludwig K.R., 2003. *ISOPLOT 3.00: A Geochronological Toolkit for Microsoft Excel*. Berkeley Geochronology Center, California, Berkeley, 39 pp.
- McCulloch M.T., Rosman K.J.R., De Laeter J.R., 1977. The isotopic and elemental abundance of ytterbium in meteorites and terrestrial samples. *Geochimica et Cosmochimica Acta*, 41(12): 1703-1707.
- Tong X.R., Liu Y.S., Hu Z.C., Chen H.H., Zhou L., Hu Q.H., Xu R., Deng L.X., Chen C.F., Yang L., Gao S., 2016. Accurate determination of Sr isotopic compositions in clinopyroxene and silicate glasses by LA-MC-ICP-MS. *Geostandards and Geoanalytical Research*, 40(1): 85-99.
- Wiedenbeck M., Alle P., Corfu F., Griffin W.L., Meier M., Oberli F., Quadt A.V., Roddick J.C., Spiegel W., 1995. Three natural zircon standards for U-Th-Pb, Lu-Hf, trace element and REE analyses. *Geostandards and Geoanalytical Research*, 19(1): 1-23.
- Woodhead J., Hergt J., Shelley M., Eggins S., Kemp R., 2004. Zircon Hf-isotope analysis with an excimer laser, depth profiling, ablation of complex geometries, and concomitant age estimation. *Chemical Geology*, 209(1-2): 121-135.
- Xu L., Hu Z., Zhang W., Yang L., Liu Y., Gao S., Luo T., Hu S., 2015. In situ Nd isotope analyses in geological materials with signal enhancement and non-linear mass dependent fractionation reduction using laser ablation MC-ICP-MS. *Journal of Analytical Atomic Spectrometry*, 30(1): 232-244.
- Zhu L., Liu Y., Ma T., Lin J., Hu Z., Wang C., 2016. In situ measurement of Os isotopic ratios in sulfides calibrated against ultra-fine particle standards using LA-MC-ICP-MS. *Journal of Analytical Atomic Spectrometry*, 31(7): 1414-1422.
- Zhu L.Y., Liu Y.S., Hu Z.C., Hu Q.H., Deng L.X., Gao S., 2013. Simultaneous analysis of major and trace elements in fused volcanic rock powders by using a hermetic vessel heater and LA-ICP-MS. *Geostandards and Geoanalytical Research*, 37(2): 207-229.
- 张春来, 刘勇胜, 高山, ZAJACZ Z., 胡兆初, 高长贵, 2011. 四合屯玄武岩斑晶中单个熔体包裹

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